



Overview

The objective of this document is to provide a quick reference for the customer to review the Masimo SafetyNet Dashboard with their team members.

The sequence of the events can be related back to any other patient monitoring system the Customer interacts with. Customers must have access, associate a patient to a device (Admit), perform actions during monitoring and dissociate the patient from the device (Discharge).




In the Masimo SafetyNet System, this is equivalent to the following:




Access

Users gain access to the Masimo SafetyNet Dashboard by having an account created for them in the “Access” tab. To create an Account the only information required is the Users email address. To create an Account, perform the following:


1. Click ‘Access’



2. Click the ‘Add User’ icon



3. Enter the Users email
 - a. **NOTE:** The field text states ‘Enter Name, email address’, however, you only enter the email address
 - b. **NOTE:** Use all lower case characters
4. Click ‘Add’


A screenshot of the "Add User" form. The form has a title "Add User" and a close button "x" in the top right corner. Below the title is a text input field containing "user@email.com" and a clear button "x" on the right. Underneath the input field are three checked checkboxes: "CareProgram Share", "Track Patients", and "Receive Notifications". At the bottom right of the form, there are two buttons: "CANCEL" and "ADD". The "ADD" button is highlighted with a green rectangular box.

Admit

Share Care Program

Patients get associated (Admitted) with the Masimo SafetyNet Portal by being sent and accepting a Care Program. To share a Care Program with a Patient, perform the following:

1. Click 'CAREPROGRAMS'




2. Locate the desired Care Program and click 'Share' on the right

COVID-19 Start



3. Enter Patient First Name
4. Enter Patient Last Name
5. Enter Patient Email or Mobile Number
6. Click 'Send'

A screenshot of the 'Share' form in the Masimo SafetyNet portal. The form is titled 'COVID-19 Start' and has two tabs: 'ONE RECIPIENT' and 'MULTIPLE RECIPIENTS'. The 'ONE RECIPIENT' tab is selected. The form contains several fields: 'First Name *', 'Last Name *', and 'Email or Mobile Number *'. The 'First Name' and 'Last Name' fields are highlighted with green boxes. The 'Email or Mobile Number' field contains the text 'patient@email.com' and is also highlighted with a green box. There are checkboxes for 'Track Progress In Dashboard' and 'Date will be set when user starts CareProgram'. An 'Optional Message for Recipient' field is also present. At the bottom right of the form, there are 'CANCEL' and 'SEND' buttons. The 'SEND' button is highlighted with a green box.


7. Patient will receive an email

Monitoring

Grid View

The Grid View provides an overview of all the patients being monitored in the one View.

- ▶ The Values displayed are the last values received from the patient
- ▶ If values have not been received in 15 minutes (default) all the parameters will be dashed out
 - ▼ This is configurable by clicking on 'Grid Setup'



- ▶ The timestamp in the top right of the box represents the last time a patient value was received
- ▶ When clicking on a patient grid, it will take you to the trends of the collected patient data

Alert Management

Patient Alerts may be accessed by clicking on the number in the 'red circle' in the top right-hand corner of the patient grid box.



This will open the list of patient alerts.

Alerts		
Messages ▲	Date ▲	Actions
SpO2 Low: 50 %	04/17/2020 04:45 pm	VIEW CLEAR


If there is any additional information configured for patients to respond with during Alert periods, you can click 'View' to get additional details. To clear the Alert from the View, click 'Clear'

List View



To access the List View, click the 'List View' icon



- ▶ All data in the List View may be sorted by clicking on any column header
 - ▼ For example, you may sort the Alerts from highest to lowest by clicking on the Alerts column
- ▶ You may access Patient Health Data (Trends) by clicking on the 'Health Data' icon




- ▶ Any Icon with the pencil icon may be edited, for example 'MRN' & 'DOB'


M.R NUMBER	D.O.B
	

- ▶ If desired you can enter notes into the 'Notes' column
 - ▼ **NOTE:** The user must click on the Notes Icon to see the Notes

Notes



View Care Program

You can view the patients Care Program, which includes their daily question responses by scrolling all the way over to the right and clicking the  icon and selecting 'View Patients CP'

Care Actions				
Onboarding ▼ ✕				
CareProgram Items ▲	Responses	Due On ▲	Last Action Date	
Please Enter Your Height	-	04/17/2020	-	
Please Enter Your Weight	-	04/17/2020	-	
Continuous Monitoring ▼ ✕				
CareProgram Items ▲	Responses	Due On ▲	Last Action Date	
✓ Oxygen Saturation (SpO2)	98 %	04/17/2020	04/17/2020 04:45:14 pm	
✓ Heart Rate	59 bpm	04/17/2020	04/17/2020 04:45:15 pm	
▼ Value within range at 04:45:17 pm ▲				

Discharging Patient

Stop Care Program

To stop a Care Program, perform the following:

1. Navigate to List View
2. Locate the patient on the List View Dashboard
3. Scroll over to the right until you can see the 'Stop Care Program' column
4. Click Stop Care Program Icon


Stop CareProgram



5. This will remove the Care Program from the patients phone the next time the Masimo SafetyNet App is synced with the Masimo SafetyNet Portal

Delete Patient from Dashboard

To remove a patient from the Dashboard, perform the following:

1. Scroll all the way over to the right and click the  icon and selecting 'Delete'

